

SMAXTEC APPLICATION EXAMPLE MASTITIS

HOW SMAXTEC DETECTS MASTITIS

Mastitis results in **considerable financial losses**, caused by **daily milk loss** in the **current and following lactation**. Depending on the time of disease occurrence, a milk loss of **2.5 kg per day** or **550 kg per lactation** can arise. **The earlier mastitis is detected, the smaller the damage.**

With smaXtec, you detect mastitis early: The **increase of body temperature** is, as for most **inflammatory diseases**, the **first indicator of a beginning disease**. According to scientific studies, body temperature increases **up to 4 days before clinical symptoms** become visible. Since the disease is detected this early, you can **act early** in



order to considerably **reduce** the number of **moderate to severe mastitis cases on your** farm.

Find below two mastitis examples: one E.coli mastitis proving how early smaXtec notifies you, and a mastitis that occurred shortly after calving which was detected early, but unfortunately didn't lead to a timely reaction.

EXAMPLE 1: E.COLI MASTITIS

Mastitis caused by E.Coli bacteria **frequently** have a **severe**, **highly acute disease progression**. E.Coli mastitis often leads to the **loss of the affected udder quarter or to the animal's death.** This is why it is **particularly important to detect E.Coli mastitis as early as possible** so you can **react** as **quickly** as possible and increase animal welfare.

Inner body temperature changes earlier than other parameters used for mastitis detection. This is the reason smaXtec is able to provide mastitis information earlier than milking systems. See for yourself in the following example: This cow had a case of E.Coli mastitis on the 12th August. smaXtec notified the farmer at 9 pm. He immediately examined the cow and took action. The milking system notified the farmer on the 13th August based on the changed conductivity of milk – more than 6 hours later!

The **severity** of disease progression in case of E.Coli mastitis leads to **considerable differences between two milkings**. In this case, there were **no signs of mastitis during the last** milking – they **only occurred afterwards**. Signs of mastitis were thus only detected during the next milking. This **delay can be fatal**, particularly when it comes to E.Coli mastitis that **requires rapid action** to **save the cow**. Once again, it shows the **high value of continuous measurement** of **inner body temperature**.





In this case, the farmer reacted **as quickly as possible** thanks to smaXtec, and was thus able to **avoid lasting damages**. The cow **recovered quickly**: Looking at the drinking cycles you can see that the cow drank considerably more water **on the 15th August**. Today, the cow is still **in the best of health**and**milk yield** was **only slightly affected** thanks to the early reaction.

But what happens if farmers do not react immediately? The following example shows how early smaXtec gives the opportunity for treatment.

EXAMPLE 2: MASTITIS AFTER CALVING

In the following example you see the **changes of the cow's body temperature** caused by a **mastitis shortly after calving**. In this case, **multiple alerts** were sent **before measures were taken**.





Here, the system sent **3 increased temperature alerts**. **Clinical signs** with flakes in the milk were only visible in the **evening of the third alert**, and this is when the **treatment** was **started**. If the cow is **carefully examined after the first notifications** by using a **California Mastitis Test**, the test already shows **slight streaks**. This leads to **earlier diagnosis** and **treatment**. In this case, the **progression** of the disease was **severe**; the **economic loss** throughout the current as well as the following lactation is **enormous**. The usage of **antibiotics** could have probably also been **avoided** if a reaction based on the smaXtec alert had taken place. Additionally, the **cow's pain** can be **relieved** earlier and thus additionally **animal welfare** is being improved.